



# Reel1.0 Quick Guide

By Frederik Holm Gjørup  
May 2021



## Contents

Welcome.....	2
Installation (Windows) .....	2
Pinning the shortcut to start/toolbar .....	2
Installation (MAC) .....	2
Getting started.....	3
Data formats .....	4
Custom .xyy file format for Reel1.0.....	4



## Welcome

Reel1.0 is a visualization tool for diffraction data and refinement results, intended for *in-situ*, *operando*, or similar large datasets. The program allows the user to “Reel” through their data and quickly assess the quality of a multitude of refinements. While the idea for the program came from parametric refinements, it is equally applicable for sequential refinements or even raw data. Reel1.0 is written in Python and is open-source and free for all to use, however, I kindly request that it is not distributed with commercial intent. I hope that you enjoy the program and encourage any and all feedback.

-Frederik

## Installation (Windows)

Reel1.0 comes as a collection of python files, and as such, require the user to install Python 3 and a few non-built-in modules. All the modules can be installed using *pip install* or *anaconda*. The latest version-requirements are: (May 2021)

Python 3.8.3

Non-built-in python modules:

Module	Version	pip install command
PyQt5	v. 5.14.2	pip install PyQt5
pyqtgraph	v. 0.11.1	pip install pyqtgraph
matplotlib	v. 3.3.3	pip install matplotlib
numpy	v. 1.19.5	pip install numpy
scipy	v. 1.6.0	pip install scipy

+ Requirements imposed by the modules. Use *pip show [module]* for information about your current version.


Once all the required modules are installed, run the *setup.bat*. This will make the *run\_Reel.bat* file and a *Reel1.0* shortcut. Open the program by double-clicking the short-cut or the *run\_Reel.bat* file. If the program does not open, check that the path to *python.exe* in *run\_Reel.bat* line 5 is valid.

Pinning the shortcut to start/toolbar

In order to pin the shortcut to the windows toolbar or start menu, right-click the *Reel1.0* shortcut, go to *properties*, add *C:\Windows\explorer.exe* in front of the *target* path, and click *OK*. Right-clicking the shortcut will now allow you to pin the shortcut to the start menu and the toolbar.

## Installation (MAC)

There are no shortcut or batch files for mac-users, instead, open a terminal window, navigate to the *Reel1.0* folder, and run the *Refinement\_evaluator\_ver1.0.py* in python.



## Getting started

New files are opened by clicking the *Open* button and selecting the desired file format in the dialog box.<sup>1</sup> Open several datasets by clicking the *Add dataset* button, and change between the opened datasets by clicking their names, the *Next* or *Previous* buttons, or with (ctrl+arrow keys). Update the current dataset by clicking *Update* or the hotkey (F5). Remove one or all dataset(s) by clicking the *Remove* buttons.

"Reel" through frames by dragging the green *Reel cursor* lines, either in the *Surface* plots or in the *Parameter* plot, or by using the arrow keys.

Select which parameters to plot and which axis to plot them on in the *Parameter plot* menu. Likewise, choose available sub plots from the *Pattern plot* menu.

You can scale the intensities separately for the *surface* plots and the *pattern* plots in the *Surface plot* and *Pattern plot* menus. Notice that scaling one will not affect the other. Adjust the colormap scaling by moving the yellow upper and lower boundaries in the *scale bar histograms* in the *surface* plot.

There are several options for zooming and panning in the plots. Change between the zoom modes for a given plot by right-clicking in the plot and selectin *Mouse Mode*. Auto-scale the range for all plots by clicking *auto-range* or using the hotkeys (A) or (space).

You can manually set the wavelength by clicking the *Set wavelength* button.

You can change the user defined default settings by editing the *ReelUserSettings.py* file in the *\_lib* folder, and restore the default settings from the *Help* menu in *Reel1.0*.



<sup>1</sup> **NB:** The file order might depend on the sorting order of the current folder!

## Data formats

Reel accepts several common data formats, but for the full range of option, use the *.xyy* format, as described below. The accepted data formats are:

- *.xyy* (custom Reel1.0 format)
- *.prf* (FullProf, prf=3)<sup>2</sup>
- *.prf* (JANA)<sup>3</sup>
- *.par* (+ associated *.fit* files, MAUD)
- *.dat* (FullProf, ins=0, 10)<sup>4</sup>
- *.xy* and *.xye* (raw data with no header)
- *.csv* (2 $\theta$  in the first line and intensities in consecutive lines)

Custom *.xyy* file format for Reel1.0

The following is a guide to the custom *.xyy* file format intended for *Reel1.0*.

The header is separated in several parts. It starts with a *mandatory* filename on the first line. Then follows a section of *keywords* ending on a colon and followed by a *value*. There are three keywords of special significance: *R<sub>wp</sub>*:, *Temperature (K)*:, and *Wavelength (Å)*:. Any additional colon-separated *keywords* and *values* are read by the program and can be plotted in the *parameter plot*.

The *keyword* section is terminated by the beginning of the comments section. The comment section is started with *COMMENTS*, and any additional lines will be ignored by the program, until the *END OF HEADER* line.

The data columns start with a column label, followed by the data. *Reel* expects the following *mandatory* columns: "*tth Y\_obs Y\_calc Y\_res Background*". Any additional columns will be added as *sub plots* with their corresponding label.

1	File name									
2	R <sub>wp</sub> : 4.702									
3	Temperature (K): 373.500									
4	Wavelength (Å): 2.52098									
5	Counts: 834837.00									
6	COMMENTS									
7	Any additional comments or metadata provided by the user									
8	END OF HEADER									
9	tth	Y_obs	Y_calc	Y_res	Background	Phase 1	Phase 2	Phase 3	Phase 4	
10	9.970	8039.000	8305.591	-266.591	8287.333	8294.666	8298.258	8287.333	8287.333	
11	10.070	7999.000	8283.805	-284.805	8264.862	8272.764	8275.903	8264.862	8264.862	
12	10.170	8045.000	8262.288	-217.288	8242.581	8251.130	8253.740	8242.581	8242.581	
13	10.270	7960.000	8241.055	-281.055	8220.490	8229.776	8231.768	8220.490	8220.490	
14	10.370	8221.000	8220.122	0.878	8198.587	8208.722	8209.986	8198.587	8198.587	
15	10.470	8110.000	8199.512	-89.512	8176.871	8187.989	8188.394	8176.871	8176.871	
16	10.570	7990.000	8179.257	-189.257	8155.342	8167.607	8166.992	8155.342	8155.342	
17	10.670	8104.000	8159.391	-55.391	8133.999	8147.613	8145.777	8133.999	8133.999	
18	10.770	7990.000	8139.962	-149.962	8112.840	8128.054	8124.748	8112.840	8112.840	
19	10.870	7894.000	8121.036	-227.036	8091.865	8108.997	8103.905	8091.865	8091.865	
20	10.970	8187.000	8102.705	84.295	8071.073	8090.531	8083.246	8071.073	8071.073	

<sup>2</sup> FullProf manual p. 73

<sup>3</sup> JANA2006 cookbook p. 773

<sup>4</sup> FullProf manual p. 74